

CHOLESTEROL

Cat. No.	Pack Name	Packaging (Content)
BLT00034	CHOL 5x50	R1: 5 x 50 ml, R2 standard: 1 x 5 ml
BLT00035	CHOL 1000	R1: 1 x 1000 ml
BLT00036	CHOL 250	R1: 1 x 250 ml, R2 standard: 1 x 5 ml

EN



INTENDED USE


Diagnostic reagent for quantitative *in vitro* determination of Cholesterol in human serum and plasma.

CLINICAL SIGNIFICANCE

Measurement of serum cholesterol levels can serve as an indicator of liver function, biliary function, intestinal absorption, propensity towards coronary artery disease, thyroid function and adrenal disease. Cholesterol levels are important in the diagnosis and classification of hyperlipoproteinaemias. Stress, age, gender, hormonal balance and pregnancy affect normal cholesterol levels.

PRINCIPLE

This reagent is based on the formulation of Allain et al and the modification of Roeschlau with further improvements to render the reagent stable in solution.



where:

CHE = Cholesterol Esterase
 CHO = Cholesterol Oxidase
 4AAP = 4-aminoantipyrine
 POD = Peroxidase

- Cholesterol esters are enzymatically hydrolysed by cholesterol esterase to cholesterol and free fatty acids.
- Free cholesterol, including that originally present, then oxidized by cholesterol oxidase to cholest-4-en-3-one and hydrogen peroxide.
- The hydrogen peroxide combines with 4-aminoantipyrine to form a chromophore (quinoneimine dye) which may be quantitated at 505 nm.

REAGENT COMPOSITION

R1

Good's Buffer 50 mmol/l
 Phenol 5 mmol/l
 4-aminoantipyrine 0.3 mmol/l
 Cholesterol esterase ≥ 200 U/l
 Cholesterol oxidase ≥ 50 U/l
 Peroxidase ≥ 3 kU/l

R2 standard See bottle label

REAGENT PREPARATION

Reagent is liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

SPECIMEN COLLECTION AND HANDLING

Use serum, plasma (heparin, EDTA).

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability

in serum / plasma: at 20–25°C 7 days
 at 4–8 °C 7 days
 at -20°C 3 months

Discard contaminated specimens.

CALIBRATION

Calibration with the standard included in the kit or calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

UNIT CONVERSION

mg/dl x 0.026 = mmol/l

EXPECTED VALUES ²

Adult

Desirable blood Cholesterol < 200 mg/dl
 Borderline high blood Cholesterol 200 – 239 mg/dl
 High blood Cholesterol > 239 mg/dl

Child

Desirable blood Cholesterol < 170 mg/dl
 Borderline high blood Cholesterol 170 – 199 mg/dl
 High blood Cholesterol > 199 mg/dl

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 4.2 mg/dl

Linearity: 695 mg/dl

Measuring range: 4.2 – 695 mg/dl

PRECISION

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	126.81	1.58	1.26
Sample 2	226.85	2.15	0.96

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	106.46	1.12	1.06
Sample 2	187.46	3.08	1.65

COMPARISON

A comparison between XL-Systems Cholesterol (y) and a commercially available test (x) using 40 samples gave following results:

$y = 0.995x - 4.59$ mg/dl

$r = 1.000$

INTERFERENCES

Following substances do not interfere:
 haemoglobin up to 5 g/l, bilirubin up to 20 mg/dl, triglycerides up to 2000 mg/dl.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent of the kit is not classified like dangerous but contains less than 0.1% sodium azide - classified as very toxic and dangerous substance for the environment.

WASTE MANAGEMENT

Please refer to local legal requirements.

ASSAY PROCEDURE

Wavelength: 500 (546) nm

Cuvette: 1 cm

	Reagent blank	Standard (Calibr.)	Sample
Reagent 1	1.00 ml	1.00 ml	1.00 ml
Sample	-	-	0.01 ml
Standard (Calibr.)	-	0.01 ml	-
Distilled water	0.01 ml	-	-

Mix and incubate 10 min. at 37 °C. Measure absorbance of the sample A_{sam} and standard A_{st} against reagent blank. The coloration is stable during one hour.

CALCULATION

$$\text{Cholesterol (mg/dl)} = \frac{\Delta A_{sam}}{\Delta A_{st}} \times C_{st}$$

C_{st} = standard (calibrator) concentration

Applications for automatic analysers are available on request.

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	End Point
Wavelength 1 (nm)	505
Wavelength 2 (nm)	670
Sample Volume (µl)	5/10
Reagent Volume (µl)	500/1000
Incubation time (min.)	5
Incubation temp. (°C)	37
Normal Low (mg/dl)	0
Normal High (mg/dl)	200
Linearity Low (mg/dl)	4.2
Linearity High (mg/dl)	695
Concentration of Standard	See bottle label
Blank with	Reagent
Absorbance limit (max.)	0.2
Units	mg/dl





Erba Lachema s.r.o., Karásek 1d, 621 00 Brno, CZ
 e-mail: diagnostics@erbalachema.com, www.erbamannheim.com


REFERENCES


1. Searcy, R.L. "Diagnostic Biochemistry" McGraw-Hill, New York, NY. 1969.
2. Tietz Textbook of Clinical Chemistry and Molecular diagnostics. Burtis, C.A., Ashwood, E. R., Bruns, D.E.; 5th edition, WB Saunders Company, 2012.
3. Flegg HM. Ann Clin Biochem. 1973; 11 : 79.
4. Richmond, W. Clin. Chem 1973; 19 : 1350-1356.
5. Allain, C.C. Poon, L.S, Chan, C.S.G, Richmond, W. and Fu, P.C. Clin Chem. 1974; 20: 470-475.
6. Roeschlau P, Bernt, E. and Gruber, W.A. Clin. Chem. Clin. Biochem. 1974; 12 : 226.
7. Henry, R. J. Clinical Chemistry: Principles and Techniques Harper & Row, Hagerstown, 1974.
8. Young, D. S. et al. Clin. Chem 1975; 21.
9. NCEP Expert Panel. Arch. Intern. Med. 1988; 148 : 36-69.
10. National Committee for Clinical Laboratory Standards. User evaluation of Precision Performance of Clinical Laboratory Devices. NCCLS; 1984, NCCLS Publication EP5-T.

SYMBOLS USED ON LABELS


 REF Catalogue Number


 Manufacturer


 See Instruction for Use

 LOT Lot Number

 CE Mark -
Device comply with
the Directive 98/79/EC


 Storage Temperature

 Expiry Date

 IVD In Vitro Diagnostics

 CONT Content

QUALITY SYSTEM CERTIFIED
ISO 9001 ISO 13485

 Erba Lachema s.r.o., Karásek 1d, 621 00 Brno, CZ
e-mail: diagnostics@erbalachema.com, www.erbamannheim.com