

# LACTATE DEHYDROGENASE-P

Cat. No.	Pack Name	Packaging (Content)
BLT00037	LDH 100	R1: 4 x 20 ml, R2: 1 x 20 ml



## INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of LDH in human serum and plasma (DGKCH method).

## CLINICAL SIGNIFICANCE

The enzyme Lactate dehydrogenase (LDH) is concentrated in heart, kidney, liver, muscle and body tissues. Consequently, damage to these results in increased serum levels of LDH. Elevated levels are associated with myocardial infarction, renal damage, hepatitis, anemias, malignancies and muscular disease or damage.

## PRINCIPLE

The LDH method is based on the recommendations of DGKCH (from pyruvate). This reagent uses pyruvate and is based on the method of Henry et al.



LDH catalyses the reduction of pyruvate to lactate oxidising reduced nicotinamide adenine dinucleotide (NADH) to NAD. The activity of LDH can be determined by the rate of decrease in absorbance at 340 nm as NAD is produced.

## REAGENT COMPOSITION

### R1

Tris Buffer (pH 7.5)	100 mmol/l
Pyruvate	2.0 mmol/l

### R2

NADH	1.66 mmol/l
------	-------------

## REAGENT PREPARATION

Reagents are liquid, ready to use.

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

### Two reagents method

Reagents are ready to use.

After opening, reagents are stable until expiry date at 2–8°C if stored at appropriate conditions, closed carefully and without any contamination.

### Monoreagent method

Mix 4 portion of Reagent R1 and 1 portion of Reagent R2.

Stability:

24 hours	at 15–25°C at dark
5 days	at 2–8°C at dark

## SPECIMEN COLLECTION AND HANDLING

Use serum, plasma (heparin, EDTA).

It is recommended to follow NCCLS procedures (or similar standardized conditions).

### Loss of activity:

within 24 hours	at 15–25°C	< 2%
within 3 days	at 2–8°C	< 8%
Stability at least 6 weeks	at –20 °C	

## CALIBRATION

Calibration with calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

## QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

## UNIT CONVERSION

U/l x 0.017 = µkat/l

## EXPECTED VALUES ?

At 37°C: 225 - 450 U/l

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 43.8 U/l

**Linearity:** 1200 U/l

**Measuring range:** 43.8 - 1200 U/l

Intra-assay precision Within run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	76.74	3.6	0.49
Sample 2	76.02	7.8	0.99

Inter-assay precision Run to run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	562.8	13.8	2.43
Sample 2	312.0	5.4	1.88

## COMPARISON

A comparison between XL-Systems LDH (y) and a commercially available test (x) using 40 samples gave following results:

$$y = 1.982x + 0.06 \text{ U/l}$$

$$r = 0.996$$

## INTERFERENCES

Following substances do not interfere:

bilirubin up to 20 mg/dl, triglycerides up to 500 mg/dl, haemoglobin up to 5.0 g/l. Significant hemolysis may increase LD concentration because of high levels of LD in the erythrocytes.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

## WASTE MANAGEMENT

Please refer to local legal requirements.

## ASSAY PROCEDURE

**Wavelength** 340 nm, Hg 334 nm, Hg 365 nm

**Cuvette** 1 cm

### Two reagents method – substrate start

	Reagent blank	Calibrator	Sample
Reagent 1	0.800 ml	0.800 ml	0.800 ml
Sample	–	–	0.020 ml
Calibrator	–	0.020 ml	–
Distilled water	0.020 ml	–	–

Mix and after 1 min. incubation (at 37°C ) add:

Reagent 2	0.200 ml	0.200 ml	0.200 ml
-----------	----------	----------	----------

Mix, incubate 1 min. at 37°C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

### Monoreagent method – sample start

	Reagent blank	Calibrator	Sample
Working reagent	1.000 ml	1.000 ml	1.000 ml
Sample	–	–	0.020 ml
Calibrator	–	0.020 ml	–
Distilled water	0.020 ml	–	–

Mix, incubate 1 min. at 37°C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

## CALCULATION

$$1. \text{LDH (U/l)} = \frac{\Delta A_{\text{sam}}/\text{min}}{\Delta A_{\text{cal}}/\text{min}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

2. Using factor:

$$\text{LDH (U/l)} = f \times \Delta A/\text{min}$$

f = factor

f = 8095 (at 340 nm)

Applications for automatic analysers are available on request.

## ASSAY PARAMETERS FOR PHOTOMETERS


Mode	Kinetic
Wavelength 1 (nm)	340
Sample Volume (µl)	10/20
Reagent Volume (µl)	500/1000
Lag Time (sec.)	60
Kinetic Interval (sec.)	60
No. of readings	3
Kinetic factor	8095
Reaction temperature (°C)	37
Reaction direction	Decreasing
Normal Low (U/l)	225
Normal High (U/l)	450
Linearity Low (U/l)	43.8
Linearity High (U/l)	2000
Absorbance Limit (Max.)	0.8
Blank with	Water
Units	U/l





#### REFERENCES


1. Searcy, R.L., Diagnostic Biochemistry, McGraw-Hill, New York, NY, 1969.
2. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
3. Henry, R.J., Chiamori N., Golub O.J., and Berkman S., Am.J. Clin. Path. 34(341), 1960.
4. Lum, G., Gambino, S.R., Am.J. Clin. Pathol. 61(108), 1974.
5. Bergmeyer, H. W., Methods of Enzymatic Analytatic Analysis, Ed.2, Verlag Chemie, 1965.
6. Young DS, Effects of Drugs on Clinical Laboratory Tests. Third Edition. 1990; 3 : 221-4.

#### SYMBOLS USED ON LABELS


 REF Catalogue Number


 Manufacturer

 See Instruction for Use

 LOT Lot Number

 CE Mark -  
Device comply with  
the Directive 98/79/EC


 Storage Temperature

 Expiry Date

 IVD In Vitro Diagnostics

 CONT Content

QUALITY SYSTEM CERTIFIED  
ISO 9001 ISO 13485

 Erba Lachema s.r.o., Karásek 1d, 621 00 Brno, CZ  
e-mail: [diagnostics@erbalachema.com](mailto:diagnostics@erbalachema.com), [www.erbamannheim.com](http://www.erbamannheim.com)