BILIRUBIN TOTAL&DIRECT

Cat. No.	Pack Name	Packaging (Content)
BLT00011	BIL T&D 200	R1: 2x50 ml (BIT), R2: 2x50 ml (BID), R3: 1x6 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Total & Direct Bilirubin in human serum and plasma.

CLINICAL SIGNIFICANCE

Bilirubin is a breakdown product of haemoglobin. Bilirubin formed in the reticulo endothelial system is transported bound by albumin to the liver. This bilirubin is water insoluble and is known as indirect or unconjugated bilirubin. In the liver, bilirubin is conjugated to glucuronic acid to form direct bilirubin. Conjugated bilirubin is excreted via the biliary system into the intestine. Here it is metabolised by bacteria to urobilinogen and stercobilinogen.

TOTAL BILIRUBIN = INDIRECT BILIRUBIN + DIRECT BILIRUBIN

Total Bilirubin is elevated in obstructive conditions of the bile duct, hepatitis, cirrhosis, in haemolytic disorders and several inherited enzyme deficiencies.

Indirect Bilirubin is elevated by pre-hepatic causes such as haemolytic disorders or liver diseases resulting in impaired entry transport or conjugation within the liver. Monitoring of indirect bilirubin in neonates is of special importance as it is the indirect (or free) bilirubin bound to albumin that is able to cross the blood brain barrier more easily increasing the danger of cerebral damage.

PRINCIPLE

Modified method of Pearlman & Lee in which a surfactant is used as a solubilizer.

Bilirubin glucuronate reacts directly with sulphodiazonium salt and forms coloured derivative azobilirubin. The colour intensity of formed azobilirubin measured at 540 - 550 nm is proportional to direct bilirubin concentration in the sample.

REAGENT COMPOSITION

R1		
Sulphanilic Acid	4.62 mmol/l	
HCl	117.6 mmol/l	
Cetrimonium Bromide	27.44 mmol/l	
R2		
Sulphanilic Acid	14.61 mmol/l	
HCl	117.6 mmol/l	
R3		
Sodium Nitrite	145 mmol/l	

REAGENT PREPARATION

Test	Volume of working Reagent	Add		
		Reagent 1	Reagent 2	Reagent 3
Bilirubin Total	10 ml	10 ml	-	0.2 ml
	25 ml	25 ml	-	0.5 ml
	50 ml	50 ml	-	1.0 ml
	100 ml	100 ml	-	2.0 ml

Test	Volume of working Reagent	Add		
		Reagent 1	Reagent 2	Reagent 3
Bilirubin Direct	10 ml	-	10 ml	0.1 ml
	25 ml	-	25 ml	0.25 ml
	50 ml	-	50 ml	0.5 ml
	100 ml	-	100 ml	1.0 ml

WORKING REAGENT STABILITY

The working Total & Direct Bilirubin reagents are stable for 7 days at 2-8°C, when protected from contamination and light. It is recommended to prepare fresh working solution before assay is performed.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2-8°C.

SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum or plasma (heparin, EDTA).

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability of total bilirubin: 1 day at 15-25°C
7 days at 2-8°C
3 month at -20°C

Stability of direct bilirubin: 2 days at 15-25°C
7 days at 2-8°C
3 month at -20°C
In case of immediate freezing.

Freeze only once!

Discard contaminated specimens.

CALIBRATION

Calibration with the factor (see ASSAY PARAMETERS table) or with calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

UNIT CONVERSION

mg/dl x 16.95 = µmol/l

EXPECTED VALUES ³

TOTAL BILIRUBIN:

Adults: 0 - 2.0 mg/dl
Cord < 2 mg/dl

Newborns, premature:

0 - 1 d 1.0 - 8.0 mg/dl
1 - 2 d 6.0 - 12.0 mg/dl
3 - 5 d 10.0 - 14.0 mg/dl

Newborns, full term:

0 - 1 d 2.0 - 6.0 mg/dl
1 - 2 d 6.0 - 10.0 mg/dl
3 - 5 d 4.0 - 8.0 mg/dl

DIRECT BILIRUBIN:

Adults and infants: 0 - 0.2 mg/dl

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Measuring range Bilirubin Total: 0.08 - 23 mg/dl

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	1.58	0.015	0.93
Sample 2	3.90	0.035	0.90

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	1.36	0.023	1.67
Sample 2	5.23	0.068	1.31

Measuring range Bilirubin Direct: 0.18 - 23 mg/dl

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	9.18	0.23	2.45
Sample 2	15.12	0.32	2.10

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	2.54	0.10	3.78
Sample 2	0.78	0.03	3.37

COMPARISON

A comparison between XL-Systems Bilirubin Total & Direct (y) and a commercially available test (x) using 40 samples gave following results:

Bilirubin Total: y = 1.091 x - 0.087 mg/dl r = 0.992
Bilirubin Direct: y = 1.029 x + 0.037 mg/dl r = 1.000

INTERFERENCES

Following substances do not interfere:

Bilirubin Total: haemoglobin up to 7.5 g/l, triglycerides up to 1500 mg/dl.

Bilirubin Direct: haemoglobin up to 10 g/l, triglycerides up to 1000 mg/dl.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. R1 contains 1% cetrimonium bromide and is classified as eye and skin irritant.

Warning

Hazard statement:

H315 Causes skin irritation.
H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.
P302+P352 IF ON SKIN: Wash with plenty of water.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

R2 contains a small amount of sulphanic acid and may cause an allergic skin reaction. R3 contains <0.1 sodium azide which is toxic and dangerous to the environment.

WASTE MANAGEMENT

Please refer to local legal requirements.

ASSAY PROCEDURE

Wavelength 546/630 (670) nm
Cuvette 1 cm

	Reagent blank	Standard (Cal.)	Sample
Working reagent	0.500 ml	0.500 ml	0.500 ml
Sample	-	-	0.025 ml
Standard (Cal.)	-	0.025 ml	-
Distilled water	0.025 ml	-	-

Mix and incubate 5 min. at 37 °C. Measure absorbance of the sample A_{sam} and standard A_{st} against reagent blank.

CALCULATION

$$1. \text{ Total \& Direct Bilirubin (mg/dl)} = \frac{\Delta A_{sam}}{\Delta A_{st}} \times C_{cal} \quad C_{cal} = \text{calibrator concentration}$$

2. Using factor:

Total & Direct Bilirubin (mg/dl) = f x ΔA
f = factor (see ASSAY PARAMETERS table)

Applications for automatic analysers are available on request.

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