




Influ-A&B Respi-Strip

Rapid diagnostic test for *in vitro* detection of Influenza A & Influenza B Viruses in swabs, washings or aspirates of nasopharyngeal secretions


- Easy to use
- Easy to read
- Specific
- Sensitive
- Reliable
- Rapid (15')
- 2x-sided sticks



www.corisbio.com

Influ-A&B Respi-Strip

Results



Performance compared to TR-FIA (n=89)

	Influenza A (side A)	Influenza B (side B)
Specificity	100% (56/56)	100% (53/53)
Sensitivity	96.9% (32/33)	97.2% (35/36)
PPV	100% (32/32)	100% (35/35)
NPV	98.2% (56/57)	98.1% (53/54)
Accuracy	98.9% (88/89)	98.9% (88/89)

For more technical information

client.care@corisbio.com


Ordering Information

Product Name	Description	Storage	Code
Influ-A&B Respi-Strip	25 tests, buffer	4°C to 30°C	C-1012
Influ-A&B Uni-Strip	10 individual strips, buffer, sample taking material (test tubes, swabs, transfert pipette, tube holder)	4°C to 30°C	K-1512

All diagnostics available in Coris BioConcept catalogue :

Adenovirus, Adenovirus (Respiratory), Adenovirus 40 and 41, *Clostridium difficile*, Bovine Coronavirus, *Cryptosporidium*, *Escherichia coli* O157, *Escherichia coli* K99, Giardia, *Helicobacter pylori*, Influenza A and B viruses, *Legionella pneumophila*, Respiratory Syncytial Virus, Rotavirus. Mefloquine, Proguanil. *Leishmania spp.*, *Trypanosoma cruzi*, *P. aeruginosa* efflux.

Distributed by:



Science Park CREALYS
Rue Jean Sonet 4A - 5032 Gembloux - Belgium
Ph.: +32 (0) 81 719 910 - Fx.: +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

MADE IN BELGIUM

Fo 7.5-32-12/05
January 2012