

# Mefloquine-Strip


**Rapid test for *in vitro* detection of Mefloquine in urine**

- Easy to read
- Rapid (10')
- Specific
- Sensitive
- No equipment required
- No dilution



[www.corisbio.com](http://www.corisbio.com)

# Mefloquine-Strip



## Performances of the test

An evaluation has been conducted in two countries on a panel of 80 individuals after intake of standard adult dose of Mefloquine for chemoprophylaxis.

Mefloquine was detected in 100% of urine specimens with the Mefloquine-Strip test.

The detection limit of Mefloquine ranges from 1 to 5 µg/mL.

## For more technical information

[client.care@corisbio.com](mailto:client.care@corisbio.com)


## Ordering Information

Product Name	Description	Storage	Code
Mefloquine-Strip	10 tests/kit	4°C to 30°C	C-10T2

### All diagnostics available in Coris BioConcept catalogue :

Adenovirus, Adenovirus (Respiratory), Adenovirus 40 and 41, *Clostridium difficile*, Bovine Coronavirus, *Cryptosporidium*, *Escherichia coli* O157, *Escherichia coli* K99, Giardia, *Helicobacter pylori*, Influenza A and B viruses, *Legionella pneumophila*, Respiratory Syncytial Virus, Rotavirus. Mefloquine, Proguanil. *Leishmania spp.*, *Trypanosoma cruzi*, *P. aeruginosa* efflux.

Distributed by:



Science Park CREALYS  
Rue Jean Sonet 4A - 5032 Gembloux - Belgium  
Ph.: +32 (0) 81 719 910 - Fx.: +32 (0)81 719 919  
[sales@corisbio.com](mailto:sales@corisbio.com) - [www.corisbio.com](http://www.corisbio.com)

MADE IN BELGIUM

Fo 7.5-32-T2/01  
January 2012