





C. diff-Strip

Rapid diagnostic test for *in vitro* detection of *Clostridium difficile* antigen in stool specimens

- Easy to read
- Easy and ready to use
- Rapid (15')
- Highly accurate
- Easy to interpret
- Specific
- Sensitive

C. diff-Strip



Performance compared to Culture (n=250)

		(95% confidence interval)
Sensitivity	98.2%	(89 to 99.9%)
Specificity	97.4%	(93.8 to 99.1%)
PPV	91.5%	(80.6 to 96.8%)
NPV	99.5%	(96.7 to 100%)
Accuracy	97.6%	

For more technical information

client.care@corisbio.com

Ordering Information

Product name	Description	Storage	Code
C. diff-Strip	25 tests, buffer	4°C to 30°C	C-1020
Clostridium K-SeT	20 cassettes, buffer	4°C to 30°C	K-1220
Clostridium K-SeT,	20 cassettes, 20 faecal sampling systems	4°C to 30°C	K-1520
with collection set			


All diagnostics available in Coris BioConcept catalogue:

Adenovirus, Adenovirus (Respiratory), Adenovirus 40 and 41, *Clostridium difficile*, Bovine Coronavirus, *Cryptosporidium*, *Escherichia coli* O157, *Escherichia coli* K99, *Giardia*, *Helicobacter pylori*, Influenza A and B viruses, *Legionella pneumophila*, Respiratory Syncytial Virus, Rotavirus.

Mefloquine, Proguanil.

Leishmania spp., *Trypanosoma cruzi*, *P. aeruginosa* efflux.

Distributed by:



Science Park CREALYS
Rue Jean Sonet 4A - 5032 Gembloux - Belgium
Ph.: +32 (0) 81 719 910 - Fx.: +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

MADE IN BELGIUM

7.5-32-20/01
September 2012