

demotek

DOUBLE BLOOD BAG 450ML

Donation No:

Collection Date:

Expiry Date:

ABO/RhD Group



DEMOPHORUS LIMITED
196 Arch. Makarios III Ave.,
2850 Limassol, Cyprus, EU,
Info@demophorus.com



	600mL
	CPDA-1 63mL
	450mL

CPDA-1 Formula:
Each 100mL contains:
Citric Acid (anhydrous) 0.299g
Sodium Citrate (dihydrate) 2.630g
Monobasic Sodium Phosphate (monohydrate) 0.225g
Dextrose (monohydrate) 3.190g
Adenine 0.0275g
Water for injection to 100mL

NEEDLE : 16G

GAUGE :

Do not vent.

For single use only.

Do not use if there is any visible evidence of deterioration.

Blood bag in unopened overwraps should be used within 15days

from opening the aluminium foil pouch and prior to the expiry date.

Sterile/non-toxic/non-pyrogenic fluid path. Steam sterilization.

Rx ONLY

UDI



2025.11

2028.10

(01)05060231201220(17)281000(10)2511303

REF BBD043



1DE0BBD043

LOT 2511303



0002511303

TRANSFER BAG

Expiry Date:



400mL

Rx ONLY



 2028.10

LOT 2511303




2DE0BBD043



0002511303

DOUBLE BLOOD BAG 450ML

4
UNITS



CPDA-1 Formula:

CPDA-1 Formula:
Each 100ml contains:

Each 100mL contains:
Citric Acid (anhydrous)

Citric Acid (anhydrous)
Monobasic Sodium Phosphate (monohydrate)

Adenine

0299a

0,299g
0,222g0,222g
0,0275g

Sodium Citrate (dihydrate)


Sodium Citrate (dihydrate)
Dextrose (monohydrate)

Dextrose (monohydrate)
Water for injection to 100ml

2.630g

2,630g
3,190g

3,190g



The plastic is of DEHP PVC. Should not be used on fetuses and neonates.

Should be used with caution on children, pregnant and lactating women.

CAUTIONS:

- CAUTIONS:**
1. Do not use if there is any visible sign of deterioration.
 2. Do not use blood bag more than 15days after opening the aluminum foil pack.



2025.11

2028.10



DEMOPHORUS LIMITED
196 Arch. Makarios III Ave.,
3030 Limassol, Cyprus, EU.
info@demophorus.com

(01)15060231201227(17)281000(10)2511303

REF BBD043



1DE0BBD043

LOT 2511303



0002511303