

TourniKidz™

Single Patient Use Tourniquet

Distraction therapy for Kidz

Combining single use benefits of infection control and employing imagery to give a less 'medical appearance', it is the perfect choice of tourniquet for paediatric healthcare professionals.



CE

Single Use Paediatric Tourniquets

tournikidz is a tourniquet that is optimised for use with children. Smaller, more comfortable, and employing imagery promoting distraction therapy, tournikidz really has put the young patient at the focus of its development.

Distraction Therapy

Colourful design, non-medical appearance draws the patient to take his or her mind away from venepuncture. The tournikidz can be given as a souvenir of their visit.

Single Patient Use

Designed for effective infection prevention.

Single handed Operation

Intuitive, safe and effective in management of blood sampling and IV cannulation.

Comfort

strong yet soft next to the skin avoiding pinching and discomfort.

tournikidz is available in a dual box containing 100 tourniquets with 50 Elly elephant and 50 Sammy snake designs.

Code: TSKKIDZ412.75DP


Instructions




Peel to expose adhesive




Pass around arm printed side up



Two hand tension



Peel off and dispose



TourniKidz
is a Trade
Mark of ASep
Healthcare Ltd

www.asephealthcare.com

T: +44 (0)1442 899942

Liverpool Science Park,
131 Mount Pleasant, Liverpool, L3 5TF

