

WHEN IT COMES TO
PATIENT COMFORT
AND SATISFACTION

think

of the ramifications



THE DALE BENDABLE ARMBORD STABILIZES AND
SECURES ARTERIAL AND PERIPHERAL LINES IN THE
OPERATING ROOM, RECOVERY ROOM, AND PATIENT
CARE UNITS.

WHEN YOU NEED SECURITY,
THINK DALE.



Dalemed.com

©2014 Dale Medical Products, Inc. All rights reserved.

DALE® BENDABLE ARMBOARD

The ArmBoard can be custom-shaped to any desired position. The ArmBoard can be shaped to achieve more prominent radial artery access. A slight bend of the ArmBoard will provide the patient with a more natural hand/wrist position. The small and x-small sizes can be shaped to accommodate an infant's wrist or heel. The ArmBoards are covered in soft, comfortable material, and fastened to the wrist/heel with adjustable straps with hook and loop closure which allows the clinician easy access to the catheter site.

- Bendable inner core
- Comfortably maintains wrist in hyperextended position
- Secures the wrist
- Improves skin care
- Improves patient comfort
- Trouble-free monitoring
- Not made with natural rubber latex.



MR Conditional, See Conditions For Use by visiting our website for detailed information on scanning safely with this device.

ORDERING INFORMATION

Product Number	Size	Quantity
H8410 650 1	Large, 23cm x 12cm (9" x 4¾")	10 per box
H8410 651 1	Medium, 14cm x 8cm (5½" x 3")	10 per box
H8410 652 1	Small, 11cm x 4cm (4½" x 1¾")	10 per box
H8410 653 1	X-Small, 11cm x 2.5cm (4¼" x 1")	10 per box



800-343-3980
www.dalemed.com

Dale® is a registered trademark of Dale Medical Products, Inc. in the USA & EU.
©2014 Dale Medical Products, Inc.

Rev D PS-650-653