


WHEN IT COMES TO  
PATIENT COMFORT  
AND SATISFACTION

# think

of the ramifications



THE DALE ABDOMINAL BINDER, SPECIFICALLY DESIGNED FOR POST-OPERATIVE USE, IS MADE ENTIRELY OF DALE'S EXCLUSIVE STRETCH MATERIAL, WHICH PERMITS HOOK AND LOOP CLOSURE AT ANY POINT FOR EVENLY DISTRIBUTED COMPRESSION. THIS STRETCH COMPRESSION ALLOWS FOR SUDDEN LUNG EXPANSION, WHICH IS NECESSARY FOR COUGHING AND DEEP BREATHING THERAPY.

**WHEN YOU NEED SECURITY,  
THINK DALE.**



[Dalemed.com](http://Dalemed.com)

©2014 Dale Medical Products, Inc. All rights reserved.

# DALE® ABDOMINAL BINDER WITH EASYGRIP™ STRIP

Studies<sup>1</sup> show that a post-operative regimen that includes the Dale Abdominal Binder with EasyGrip Strip encourages better pulmonary function, allows earlier ambulation, and alleviates incisional pain.

- EasyGrip Strip holds drainage bulbs in place.
- Evenly distributes compression.
- Facilitates deep breathing.
- Not made with natural rubber latex.

## THE DALE BINDER IN THE POST-OPERATIVE REGIMEN

The Dale Abdominal Binder with EasyGrip Strip splints the wound and encourages patients to breathe deeply, turn in bed, cough to clear their lungs, and ambulate immediately following surgery. As a result, there is less risk of sedentary complications such as atelectasis, pneumonia, and phlebitis.

## APPLICATIONS

The Dale Abdominal Binder with EasyGrip Strip provides support and comfort following procedures such as:

- Cesarean Section
- Gastroplasty
- Hysterectomy
- Abdominoplasty
- Liposuction
- Childbirth



HCPCS Code A4649

<sup>1</sup>Kuan TA, Serrett C, and Ali J; The Effect of Abdominal Binder of Postoperative Pulmonary Function. *Infections in Surgery* 1983;2(11):875-883

## ORDERING INFORMATION

Product Number	Size	Quantity
H84104101	9" Wide, Stretches to fit 76-114cm (30-45")	1 per box
H84104111	9" Wide, Stretches to fit 117-157cm (46-62")	1 per box
H84104181	9" Wide, Stretches to fit 152-191cm (60-75")	1 per box
H84108101	12" Wide, Stretches to fit 76-114cm (30-45")	1 per box
H84108111	12" Wide, Stretches to fit 117-157cm (46-62")	1 per box
H84108181	12" Wide, Stretches to fit 152-191cm (60-75")	1 per box
H84108201	12" Wide, Stretches to fit 183-239cm (72-94")	1 per box
H84109201	15" Wide, Stretches to fit 183-239cm (72-94")	1 per box



800-343-3980  
www.dalemed.com

Dale® is a registered trademark of Dale Medical Products, Inc. in the USA & EU.  
EasyGrip™ is a trademark of Dale Medical Products, Inc. in the USA.  
©2014 Dale Medical Products, Inc.