

WHEN IT COMES
TO SECURING
NASOGASTRIC TUBES

think

of the ramifications



THE DALE NASOGASTRIC TUBE HOLDERS ARE DESIGNED TO REDUCE PATIENT COMPLICATIONS AND THE COST AND TIME ASSOCIATED WITH TUBE DISLODGEEMENT.

THE UNIQUE DUAL INTER-LOCKING TABS SPIRAL AROUND THE NASOGASTRIC TUBE ENSURING SUPERIOR TUBE SECURITY. THE NON-ADHERENT PLASTIC TIPS ON BOTH TABS MAKE TUBE REMOVAL SAFER AND EASIER WITHOUT SCISSORS.

WHEN YOU NEED SECURITY,
THINK DALE.



Dalemed.com

©2014 Dale Medical Products, Inc. All rights reserved.

DALE® NASOGASTRIC TUBE HOLDERS

The Dale NasoGastric Tube Holders are ideal for holding nasogastric feeding and aspiration tubes such as: Salem Sump® tubes, Levin stomach tubes, nasoenteric feeding tubes and naso-gastric intestinal tubes.

The holder is an improvement over tape and other types of adhesive backed holders, because it features a skin-friendly adhesive that adheres to the skin of the nose quickly and comfortably. The holder can remain in place for up to 72 hours. The holder centers the tube exiting either from the left or right nostril, allowing the tube to “hang from” of the nose in a more natural position.

- Easy to use.
- Securely holds all nasogastric tubes in place.
- Prevents nasal irritation, erosion or necrosis.
- Helps prevent inadvertent dislocation or removal of tube.
- Provides maximum tube security.
- Nose pad stretches and conforms to various nose shapes.
- Non-adherent tips allow holder to be removed easily without sticking to exam gloves.
- Not made with natural rubber latex.



MedAssets Contract
Number MS02828

ORDERING INFORMATION

Product Number	Size	Quantity
H84101601	Dale® NasoGastric Tube Holder One size fits most	50 per box
H84101611	Dale® NasoGastric Tube Holder Small	50 per box



800-343-3980
www.dalemed.com

Dale® is a registered trademark of Dale Medical Products, Inc. in the USA & EU.
Salem Sump® is a registered trademark of Covidien AG Corporation in Switzerland.
©2014 Dale Medical Products, Inc.