


WHEN IT COMES
TO SECURING
TRACHEOSTOMY TUBES

think

of the ramifications



DALE TRACHEOSTOMY TUBE HOLDERS ARE DESIGNED TO PROVIDE SECURE POSITIONING AND TO MINIMIZE MOVEMENT OF THE TRACHEOSTOMY TUBE FOR ALL YOUR PATIENTS. WHETHER YOU CHOOSE THE WELL-KNOWN DALE 240 BLUE® FOR ADULTS, PEDIPRINTS® SERIES FOR KIDS, OR THE BARIATRIC 244 BLUE™ YOU CAN BE SURE IT'S THE RIGHT CHOICE.



Dalemed.com

©2014 Dale Medical Products, Inc. All rights reserved.

DALE® TRACHEOSTOMY TUBE HOLDERS

240 BLUE®, 244 BLUE™ & PEDIPRINTS®

Unnecessary movement of the trach tube can result in accidental dislodgement or displacement of the tube, tracheoesophageal fistula, tracheal stenosis, or airway granuloma. That's why it is important to know that no other securement method stabilizes tracheostomy tubes better than Dale®. Now you can eliminate the frustrations associated with twill ties and other holders while minimizing secondary complications at the same time.

- Secure fastener tabs limit movement of the tracheostomy tube reducing accidental decannulation.
- Exclusive moisture repellent neckband lining reduces the risk of skin breakdown by minimizing the amount of moisture absorbed.
- Adjustable two piece neckband.
- PediPrints® series offers kids fun, comforting prints.†
- Fastener tabs quickly and easily attach to any size tracheostomy tube.
- Stretch material‡ allows for cough reflex, accommodates edema, and ensures a snug fit without compromising security.
- Not made with natural rubber latex.

† On the 241 and 242 tube holder

‡ On the 240 Blue® tube holder



HCPCS Code A7526

ORDERING INFORMATION

Product Number	Size	Quantity
H8410 240 1	240 Blue® 1" wide band, fits up to 19.5" neck	10 per box
H8410 241 1	241 PediPrints® 1" wide band, fits up to 18" neck	10 per box
H8410 242 1	242 PediPrints® ¾" wide band, fits up to 9" neck	10 per box
H8410 244 1	244 Blue™ 1½" wide band, fits up to 28" neck	10 per box



800-343-3980
www.dalemed.com

Dale®, 240 Blue®, PediPrints® are registered trademarks of Dale Medical Products, Inc. in the USA & EU.

244 Blue™ is a trademark of Dale Medical Products, Inc. in the USA.

©2014 Dale Medical Products, Inc.