

SINGLE BLOOD BAG 450ML

Donation No:

Collection Date:

Expiry Date:

ABO/RhD Group



	600mL
	CPDA-1 63mL
	450mL

CPDA-1 Formulas:	
Each 100mL contains:	
Citric Acid (anhydrous)	0.299g
Sodium Citrate (dihydrate)	2.630g
Monobasic Sodium Phosphate (monohydrate)	0.222g
Dextrose (monohydrate)	3.190g
Adenine	0.0275g
Water for injection to 100ml	

NEEDLE : 16G

Do not vent.

Do not vent.
For single use only.

Do not use if there is any visible evidence of deterioration.

Blood bag in unopened overwraps should be used within 15 days from opening the aluminium foil pouch and prior to the expiry date. Sterile, non-toxic, non-pyrogenic fluid path. Steam sterilization.

Rx ONLY



2025.11

2028.10

(01)05060231200193(17)281000(10)2511302

REF BBS040



1DE0BBS040

LOT 2511302




0002511302

demotek

SINGLE BLOOD BAG 450ML

5
UNITS




NEEDLE • 16G
GAUGE

CPDA-1 Formula:

Each 100mL contains:

Citric Acid (anhydrous)	0,299g	Sodium Citrate (dihydrate)	2,630g
Monobasic Sodium Phosphate (monohydrate)	0,222g	Dextrose (monohydrate)	3,190g
Adenine	0,0275g	Water for injection to 100ml	



The plastic is of DEHP PVC. Should not be used on fetuses and neonates. Should be used with caution on children, pregnant and lactating women.

CAUTIONS:

1. Do not use if there is any visible sign of deterioration.
2. Do not use blood bag more than 15days after opening the aluminum foil pack.

UDI



2025.11
2028.10

DEMOPHORUS LIMITED
196 Arch. Makarios III Ave.,
3030 Limassol, Cyprus, EU,
info@demophorius.com

(01)15060231200190(17)281000(10)2511302

REF BBS040

LOT 2511302



1DE0BBS040



0002511302