

demotek

TRIPLE BLOOD BAG 450ML

Donation No:


Collection Date:

Expiry Date:

ABO/RhD Group



DEMOPHORUS LIMITED
196 Arch, Makarios III Ave.,
2050 Limassol, Cyprus, EU.
info@demophorus.com



CPDA-1 Formula:
Each 100mL contains:
Citric Acid (anhydrous) 0.299g
Sodium Citrate (dihydrate) 2.630g
Monobasic Sodium Phosphate (monohydrate) 0.225g
Dextrose (monohydrate) 3.190g
Adenine 0.0275g
Water for injection to 100ml

NEEDLE : 16G

GAUGE :

Do not vent.

For single use only.

Do not use if there is any visible evidence of deterioration.

Blood bag in unopened overwraps should be used within 15 days

from opening the aluminium foil pouch and prior to the expiry date.

Sterile, non-toxic, non-pyrogenic fluid path. Steam sterilization.

Rx ONLY

UDI



2025.11

2028.10

(01)05060231200209(17)281000(10)2511304

REF BBT040



1DE0BBT040

1

LOT 2511304



0002511304

demotek

TRANSFER BAG

Donation No:

Collection Date:

Expiry Date:

ABO/RhD Group



400mL

5

store platelets up to 5 days at 22°C ± 2°C

Do not vent.
For single use only.
Do not use if there is any visible evidence of deterioration.
Sterile, non-toxic, non-pyrogenic fluid path. Steam sterilization.

Rx ONLY

2025.11
2028.10



DEMOPHORUS LIMITED
196 Atch, Makarios III Ave.,
2030 Limassol, Cyprus. EU.
info@demophorus.com

REF BBT040



2DE0BBT040

LOT 2511304



0002511304

2

demotek


TRANSFER BAG

Donation No:

Collection Date:

Expiry Date:

ABO/RhD Group



Do not vent.
For single use only.
Do not use if there is any visible evidence of deterioration.
Sterile, non-toxic, non-pyrogenic fluid path. Steam sterilization.

Rx ONLY



DEMOPHORUS LIMITED
196 Arch, Makarios III Ave.,
2803 Limassol, Cyprus, EU
info@demophorus.com

2025.11

2028.10

REF BBT040



3DE0BBT040

3

LOT 2511304




0002511304

TRIPLE BLOOD BAG 450ML

The diagram illustrates the experimental setup for the three-chambered system. A hand is shown adding a substance to a reservoir. The reservoir is connected to three chambers. The first chamber contains 600 mL of CPDA-1 and 63 mL of a substance. The second chamber contains 400 mL of a substance. The third chamber contains 400 mL of a substance.

CPDA-1 Formula:			
Each 100mL contains:			
Citric Acid (anhydrous)	0.299g	Sodium Citrate (dihydrate)	2.630g
Monobasic Sodium Phosphate (monohydrate)	0.222g	Dextrose (monohydrate)	3.190g
Adenine	0.0275g	Water for injection to 100ml	



1. Do not use if there is any visible sign of deterioration.
2. Do not use blood bag more than 15days after opening the aluminum foil pack.



 DEMOPHORIOUS LIMITED
196 Arch. Makarios III Ave.,
3030 Limassol, Cyprus, EU.
info@demophorius.com

REF BBT040

LOT 2511304



1DE0BBT040



0002511304