



Pressure Infusion Infu-Surg®

The latex free Infu-Surg® Pressure Infusion Bag is a safe, high-quality, reliable device for intravenous and intra-arterial infusion treatment, including A-line pressure monitoring. Designed for single patient use, the Infu-Surg may help prevent and control cross-contamination and hospital acquired infections. The Infu-Surg is a cost-effective, trusted product that medical professionals have relied on for decades. The dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole.

For detailed product information, visit www.ethoxmedical.com or contact your Ethox Medical Sales Representative.



Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use

Infu-Surg® Pressure Infusion Bag

- Pressure relief valve prevents over inflation (330 mmHg pressure relief)
- Large, oval shaped bulb allows for quick and easy inflation of bladder
- Single-handed inflation and deflation design makes it easy to use and requires minimal training
- Dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole.
- Suitable for use with external inflation sources
- Color-coded gauge makes for accurate pressure monitoring (0-300 mmHg)
- Three-way stopcock ensures precise control of pressure
- Incredibly reliable – 100% tested
- Loads quickly and easily
- Latex free, single patient, non-sterile

Cat#	Size	Description	Box
4005H	500cc	Infu-Surg with Hook	5
4010H	1000cc	Infu-Surg with Hook	5
4030	3000cc	Infu-Surg	5



Pressure Infusion Infu-Surg® Clarity

The latex free Infu-Surg® Clarity Pressure Infusion Bag is an innovative, user-friendly device for intravenous and intra-arterial infusion treatment. Designed with a clear panel, the patient label and barcode on the solution bag can be easily scanned. The dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole.

For detailed product information, visit www.ethoxmedical.com or contact your Ethox Medical Sales Representative.



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use



Welded instead of sewn for added strength and durability

360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation



Clear panel enables full visibility of solution bag to monitor fluid level and easily scan patient barcode

Infu-Surg® Clarity Pressure Infusion Bag

- Transparent panel for clear visibility of solution bag
- Dedicated hook enables quick change of solution bag while providing proper fluid bag positioning in pocket
- Pressure relief valve prevents over-inflation (330 mmHg)
- Color-coded gauge provides accurate pressure monitoring (0-300 mmHg)
- Single-handed inflation and deflation with oval-shaped bulb makes it easy to use and requires minimal training
- Welded panel and hook construction are more durable than sewn designs to withstand repeated use
- Suitable for use with external inflation sources
- Three-way stopcock ensures precise control of pressure
- Incredibly reliable – 100% tested
- Latex free, single patient, non-sterile

Cat#	Size	Description	Box
4050	500cc	Infu-Surg Clarity	5
4100	1000cc	Infu-Surg Clarity	5