

EYEGARD™

EyeGard™ is the safe, effective way to protect the eyes during surgery.



EyeGard™ is a transparent barrier to cover the closed eye during anesthesia. EyeGard™ replaces tape, eliminating difficult tape handling with a gloved hand.

EyeGard™ is a new, simple, unique and patented product consisting of a non-sterile, non-latex hypoallergenic material combined with a breathable, waterproof adhesive seal. The device is oval in shape, flexible to fit the patient, and has a tab (without adhesive) for easy application and removal.

- **Cost effective**
- **Reduces lash removals**
- **Minimizes disease transmission**



Pediatric size available

Latex Free

Order#	Description	Qty
S2020	Adult	50 pairs per box
<i>New!</i> S2020-P	Pediatric	50 pairs per box
<i>New!</i> S2020-M	Adult less sticky	50 pairs per box

SHARN INC.
ANESTHESIA

4517 George Road • Suite 200 • Tampa FL. 33634
813-889-9614 • 1-800-325-3671 • FAX: 813-886-2701
www.sharn.com • email: mailbox@sharn.com