

Genrui

EXCEED IN QUALITY AND MORE

FA. 50

Quantitative Immunoassay Analyzer





50

Quantitative Immunoassay Analyzer

Easy to use

Display: 5.0-inch touch screen
Calibration mode: Magnetic card
Maintenance-free

Rapid detection

Fast mode and standard mode
Rapid detection: Results available within 3-15 min
Measurement time: 10 sec


Technical specifications

Principle: Immunofluorescence
Sample type: Whole blood, serum, plasma
Sample volume: 100ul
Storage: 50,000 results storage
Output: Internal thermal printer
Support LIS and barcode scanner
Update: 2 USB ports
Work temperature: 10~30°C
Power supply: 100~240 AC, 50/60Hz
Net weight: 4Kg

Test procedure



Magnetic card calibration



Pipette sample



Insert cassette




Result

Test menu

Hormone	β -HCG PRL LH FSH *Prog *Testo *AMH
Thyroid Function	TSH T3 T4 *FT3 *FT4
Diabetes	HbA1c
Inflammation	PCT CRP+Hs-CRP IL6
Cardiac Marker	cTnI Myo CK-MB D-Dimer NT-proBNP
Others	FER *VitD *PSA *f-PSA *AFP *CEA

More items coming soon ...



◆ Hot-selling testing
**SARS-CoV-2 Neutralizing Antibody
Test Kit (Immunofluorescence)**

All rights reserved. Specifications subject to change without prior notice

Version 2.1



Genrui Biotech Inc.
Add: 4-10F, Building 3, Geya Technology Park, Guangming District, 518106, Shenzhen, China
Tel: +86 755 26835560/1/2 Fax: +86 755 26678789
Website: www.genrui-bio.com