


## AIRWAY MANAGEMENT

SunMed

## INTRODUCER/BOUGIE

LATEX FREE | STERILE | SINGLE USE | DISPOSABLE | CE MARK

SunMed's latex free introducers are both firm and flexible for ease of insertion. They are calibrated for accurate entry depth and available in adult and pediatric sizes.



- Low density polyethylene provides proper stiffness for ease of insertion
- Calibrated to ensure accurate distance insertion
- Pediatric size features calibrated markings on both ends allowing it to be used as a straight or coudé tip
- Pediatric size fits in 4mm to 6mm tubes
- Adult size fits in 6mm to 11mm tubes

## ENDOTRACHEAL TUBE INTRODUCER

CAT#	SIZE	DESCRIPTION	PK
9-0211-70	10Fr x 70cm	Pediatric, Coudé Tip	10
9-0212-70	15Fr x 70cm	Adult, Coudé Tip	10
9-0212-72	15Fr x 70cm	Adult, Straight Tip	10



## MALLEABLE ENDOTRACHEAL TUBE INTRODUCER

- Combines benefits of an introducer and stylette into one unit
- Pediatric size features calibrated markings on both ends allowing it to be used as a straight or coudé tip

CAT#	SIZE	DESCRIPTION	PK
9-0213-82	10FR x 70cm	Pediatric Malleable, Coudé & Straight Tip	10
9-0213-92	15FR x 70cm	Adult, Malleable, Coudé Tip	10




AIRWAY MANAGEMENT

SunMed

BOUGIE-TO-GO

LATEX FREE | STERILE | SINGLE USE | DISPOSABLE | CE MARK

SunMed's latex free Bougie-To-Go™ introducer was developed with first responders in mind. Designed with the same firm and flexible quality as our traditional introducer, the Bougie-To-Go features compact packaging. It fits in paramedic bags, pockets and smaller storage spaces.



BOUGIE-TO-GO™ ENODOTRACHEAL INTRODUCER

- Innovative compact packaging designed specifically for use by EMS professionals
- Low density polyethylene provides proper stiffness for ease of insertion
- Calibrated to ensure accurate distance insertion
- Adult size fits in 6mm to 11mm tubes

CAT#	SIZE	DESCRIPTION	PK
9-0212-82	15FR x 60cm	Adult, Coudé Tip	10



Unwraps without distorting introducer shape; keeps tip in correct insertion position

