



Humasis

Handy U-AQ

Semi-Automated Single-test Analyzer

R-50S

- Light/Portable analyzer
- Colorful 4.3 inch touch LCD display
- Auto-calibration
- Quick result at only 60 Sec.
- Possibility with 8 x battery 1.5V
- USB interface
- Connectivity to LIS(Laboratory Information System)

Technical Specification

Throughput	45 tests / hour
Memory	20,000 tests
Net weight	0.4 kg
Optional accessories	External printer, Bluetooth



Product Features

- Simple operation by Touch LCD
- Light weight and compact (portable) size
- Automatic display of test results and relevant information on LCD screen
- Result can be obtained after 60 sec.
- Compatible with various kinds of CYBOW strips (incl. Microalbumin Test)



Technical Specification

Type : Reflectance Photometer

Throughput : 45 Tests (Max.120 Tests) / Hour

Interface : USB

Dim. / Net W'gt : 188(W)x74(D)x77(H) mm / 0.40 kg

Power : 12V DC/ 3.33A or battery 1.5V x 8EA

Memory: 2000 Tests

ID input : via 4.3 inch Touch LCD

Software : User PC s/w is available

Printing : via PC s/w or External Thermal Printer

Compatible strips : CYBOW 10, 11, 2AC, 10AC etc.

Basic Accessories

AC Adaptor, Power Cable, Operation Manual, USB cable, User PC s/w

Optional Accessories

External Printer, Bluetooth

Humasis

Manufacturer : Humasis

Rm. 504 SHINWON VISION TOWER, 88 Jeonpa-ro,
Dongan-gu, Anyang-si, Gyeonggi-do, 14042, KOREA
Tel. +82-31-478-8590 / Fax, +82-31-478-8586

www.humasis.com