



THINK

DALE®

Dale® UroDynamic Tube Holder

The **Dale UroDynamic Catheter Tube Holder** is ideal for holding urodynamic catheters in position during the urodynamic monitoring process to help you achieve an accurate reading.

Key Features

- Skin friendly adhesive
- Dual inter-locking tabs
- Non-adherent adhesive tips

Catheter Fixation

Secure catheter fixation is essential in maintaining good transducer signal quality and accuracy of testing. Urodynamic catheter can easily slide out of position if not securely fixated because there is no internal balloon retention.

Current methods of using standard hospital tape are inconvenient, tend to cause skin irritation, and may compress the urethra. The Dale Urodynamic Holder features two locking tabs that secure the urodynamic catheter in place. It is custom made with skin friendly adhesive and is latex free. The Dale Urodynamic Catheter can save clinicians time, allows for precise catheter positioning and secure catheter fixation which are all essential to maintain signal quality and test accuracy. The product is ideal for all male patients.

Ordering Information

Product #	Description
#H84104901	Dale UroDynamic Tube Holder One size fits all

*Individually packaged. For single patient use.
Sold in boxes of 50*

Latex-Free. In response to clinical concerns regarding patient hypersensitivity, Dale products contain no latex.

Dale products are available from your preferred hospital/medical supply dealer.

800-343-3980 • www.thinkdale.com

Dale® is a registered trademark of Dale Medical Products, Inc. in the USA & EU.
©2011 Dale Medical Products, Inc.

Rev A ps-490



EDUCATION
INNOVATION
PERFORMANCE
RELIABILITY
VALUE