

PATIENT COMFORT AND
SATISFACTION SHOULD
NEVER BE COMPROMISED

think

of the ramifications

DALE® ABDOMINAL BINDER WITH EASYGRIP™ STRIP

CLINICAL ADVANTAGES

Features	Benefits	Clinical Advantages
Improves Patient Comfort	<ul style="list-style-type: none"> Replaces need to splint the abdomen with a pillow or other methods Replaces pillow and other methods to splint radical abdominal surgeries such as abdominoplasty, liposuction, etc. Hook and loop closures allow for easy adjustments With this added security the patient is able to get out of bed quicker, thus increased healing time, and quicker out of the hospital. Saves hospital money by decreasing length of hospital stay 	<ul style="list-style-type: none"> Encourages increased activity to decrease post-op complications such as atelectasis and pneumonia Encourages patient compliance with post-op instructions to cough, deep breathe, and ambulate Promotes recovery by increasing patient confidence and sense of security
Exclusive Stretch Material	<ul style="list-style-type: none"> Hooks can be engaged at any point along the binder for desired compression and perfect fit Can be used as an abdominal dressing holder Permits quick and easy access to wound site without removing and reapplying tape Faster to change than "Montgomery Straps" or tape Will not fray or unravel if cut to allow for tubing 	<ul style="list-style-type: none"> All-elastic material allows total circumferential compression Allows clinician to control amount of pressure applied to the abdomen Supports position of taping to enhance adherence of dermal layer Can be cut to accommodate drainage tubes Allows for quick and easy wound assessment Eliminates use of tape and associated trauma from tape removal
EasyGrip™ Strip	<ul style="list-style-type: none"> Holds up to four 100cc drainage bulbs Drainage Bulbs can be emptied or removed quickly without the hassle of removing pins or tape Allows for easier monitoring of drainage bulb contents Quick tape-free application Trouble-Free monitoring Allows for easier ambulation when drainage bulbs are in place 	<ul style="list-style-type: none"> One strip can accommodate up to four bulbs No improvised securement method Clinician can easily see contents of bulb Encourages ambulation Saves nurses time
Paneled Construction	<ul style="list-style-type: none"> Only stretches in one direction preventing riding, roping and folding over 	<ul style="list-style-type: none"> Allows full range of motion Provides firm support
Innovative Joining Stitch	<ul style="list-style-type: none"> Gap-Free stitching prevents skin from migrating through stitching Will not pinch the skin 	<ul style="list-style-type: none"> Reduces the risk of abrasion to the surgical site
Convenient Sizing	<ul style="list-style-type: none"> Three widths (9", 12" and 15") are available in lengths from 30" to 94" The #411 and #811 models can be used as one-size-fits all binder, to further reduce inventory 	<ul style="list-style-type: none"> Assures correct compression and fit Developed to accommodate bariatric surgical needs
Not Made With Natural Rubber Latex	<ul style="list-style-type: none"> First Latex-Free Abdominal Binder on the market Eliminates patient and care-giver latex reactions 	<ul style="list-style-type: none"> Prevents latex hypersensitivity



800-343-3980
www.dalemed.com

Dale® is a registered trademark of Dale Medical Products, Inc. in the USA & EU. EasyGrip™ is a trademark of Dale Medical Products, Inc. in the USA. ©2014 Dale Medical Products, Inc. CA-410_Rev_D