

WHEN IT COMES TO
PATIENT COMFORT
AND SATISFACTION

think

of the ramifications



THE DALE TRANSDUCER HOLDER PROVIDES A LEVEL
AND SECURE PLATFORM, WHICH CAN ACCOMMODATE
UP TO THREE TRANSDUCERS FOR YOUR PATIENTS ON
HEMODYNAMIC MONITORING.

WHEN YOU NEED SECURITY,
THINK DALE.



Dalemed.com

©2014 Dale Medical Products, Inc. All rights reserved.

DALE® TRANSDUCER HOLDER

The stability of the Dale Transducer Holder will support correct functioning of corresponding numerical display when positioned correctly at phlebostatic axis. The Dale Transducer Holder reduces the likelihood of inaccurate measurements due to improper positioning or leveling of transducers.

The platform provides the base in which the three level holding tabs are set. The platform is made of a soft, comfortable material, and secured to the arm with an elasticized strap with adjustable hook and loop closure.

The three holding tabs reduce the time and trouble of manually leveling and taping individual transducers to the patient's chest or arm or mounting on the bedside pole. The difficulties in taping the transducers to the patient or repositioning the bedside pole because of patient movement are eliminated with the application of the Dale Transducer Holder.

- Stabilizes up to three transducers
- Improves patient monitoring
- Quick and easy to apply and remove
- Saves nursing time
- One size fits most
- Not made with natural rubber latex.



US Patent #7,284,729 & 7,284,730

ORDERING INFORMATION

Product Number	Size	Quantity
H8410 6201	Stabilizes up to three transducers	10 per box



800-343-3980
www.dalemed.com

Dale® is a registered trademark of Dale Medical Products, Inc. in the USA & EU.
©2014 Dale Medical Products, Inc.
Rev D PS-620

OTHER DALE PRODUCTS OF INTEREST

Dale® Bendable ArmBoards

H8410 6501	Large, 23cm x 12cm (9" x 4¾")	10 per box
H8410 6511	Medium, 14cm x 8cm (5½" x 3")	10 per box
H8410 6521	Small, 11cm x 4cm (4½" x 1¾")	10 per box
H8410 6531	X-Small, 11cm x 2.5cm (4¼" x 1")	10 per box

